




Seltis

SPECIALIST IN MEDICAL LAB FREEZERS



A perfect partner for your multi cold storage purposes in a contemporary lab set up. Microprocessor based temperature controls provide outstanding temperature uniformity with high quality compressor for energy efficiency. Safety features and durable interiors provide you with the most dependable performance. Used in Laboratories, Routine Cold Storage, Biologicals, Pharmaceuticals, Clinical, Research.

- Micro processor control LED digital display with adjustable temperature range from -10°C to -25°C (0.1°C increment)
- CFC-Free polyurethane foam insulation.
- Upright type aluminum spray plate interior and painted steel board exterior.
- Audio and visual alarm.
- Temperature measuring holes are deployed on the side for user convenience.
- The built in door gasket.
- The door is equipped with a lock system for sample safety.
- Caster design for easy movement.
- Drawers are available for sample store.
- Chemically suitable CFC refrigerant are used.
- ISO 9001, CE Certified, FDA, UL, IEC 60601-1 & IEC 61010-1, ROHS



SPECIFICATIONS

Model	ST-25V410
Capacity	410L
Temperature Range	-10°C to -25°C
Ambient Temperature	10°C~+32°C
Defrosting	Manual
Cooling Type	Direct Cooling
Refrigerant	CFC FREE
Internal Dimension	645X528X1641mm
Overall Dimension	775x720x1900mm
Package Size	790x730x1950mm
Drawers/Shelves	6/6
Controller	Microprocessor, Digital display
Sensor	PT-100
Compressor Brand/QTY	EMBRACO/1
Noise Level	46 db
Alarm	High/low temperature, Sensor malfunction, Temperature controller malfunction, Door ajar
Power Consumption (Kwh/24hrs)	1.5
Power Supply	220V, 50Hz

Manufactured by

SELECT REFRIGERATION

E-283, Patel Garden, Dwarka Mor, New Delhi - 110059 INDIA

Ph : 91-11-65642376, 9717569994, 9717569992 Fax : 91-11-43852133

E-mail : info@selectrefrigeration.in

Website : www.seltis.in